

## Irrigation

### Irrigating Cannula

Luer hub, angled  
10-2077/01 27 gauge




10-2077/02 30 gauge



### Irrigating/Aspirating Cannula

Luer hub, straight  
side-by-side front openings  
includes flexible tubing and hub  
10-2078




### Aspiration Handle

hollow center, one end attaches to plain  
surgical tubing while  
the other end is  
tapered to mate with Luer needles.  
10-2079



### Silicone Irrigator

The flat base allows the irrigator to be placed  
upright on the surgical tray when not in use to avoid  
spillage and loss of fluid.  
with Luer adapter  
10-5452




### Irrigating Cannula

An olive tip helps reduce the chance of accidental  
damage to delicate tissue structures during probing  
and irrigating.  
with 1mm diameter olive tip  
10-5450 25 gauge




### Hydrocellulose Sponge

for absorbing small amounts of fluids at the  
operative site  
sterile packed, 20 per package  
10-2084/01 triangular shaped  
10-2084/02 rectangular shaped



### Sponge Holder

18cm (7")  
10-2084/03



# MICROSURGERY